



## Flame Photometer TRFP-603

[www.chinascientific.com](http://www.chinascientific.com) | [info@chinascientific.com](mailto:info@chinascientific.com)

## Description

Flame Photometer TRFP-603 is equipped with a 7-inch LCD touch screen and having with a programmed calculation of the correlation coefficient. It has the ability of simultaneous measurement of the elements Na, K, Ca, Ba, and Li. It is a flameout protection device and comes with an air compressor. It has a compassion range of elements. In addition, the repeatability of elements is 1% CV for 20 consecutive samples and can accept the sample 6L/min.

## Features:

It has a 7-inch LCD touch-screen

Consistent, simultaneous measurement of elements

It can automatically calculate the coefficient correlation

Flameout protection device

LPG is a flammable gas

Information of direct concentration

Flame size pre-selection

Range shifts were quantified

The device contains an air compressor

## Specifications :

<b>Gross Weight</b>	17 kg
<b>Linearity</b>	0.195ppm?0.39 to 3.12?ppm, 0.69ppm?1.15 to 9.2?ppm, 0.15ppm?0.25 to 5?ppm
<b>Package</b>	570×530×380mm 0.115M3
<b>Range</b>	0 to 100 ppm
<b>Repeatability</b>	1 % CV for 20 consecutive samples
<b>Respond time</b>	8s
<b>Sample uptake</b>	6L/min
<b>Sensitivity</b>	0.1ppm

[www.chinascientific.com](http://www.chinascientific.com)

W Parkway St Denton TX 76201, USA

Email: [info@chinascientific.com](mailto:info@chinascientific.com) | Website: [www.chinascientific.com](http://www.chinascientific.com)