



## Digital Heating Water Bath TRHB-603

[www.chinascientific.com](http://www.chinascientific.com) | [info@chinascientific.com](mailto:info@chinascientific.com)

## **Description**

Digital Heating Water Bath TRHB-603 is a microprocessor controller with an inner digital display showing a temperature range of (RT+5) -300 °. It holds with a capacity of 10 L. It is designed with a stainless steel interior with an aluminum ring cover.

## **Features:**

It is equipped with a microprocessor control

mechanism with a stainless steel interior

Designed with a digital display inner temperature

It is covered with an aluminum ring

It is a thermostatic device

## Specifications :

<b>Capacity</b>	10 L
<b>Dimension</b>	37x37x24 cm
<b>Power</b>	1000 W
<b>Temperature Fluctuation</b>	±1 ?
<b>Temperature Homogeneity</b>	±1 ?
<b>Temperature Range</b>	(RT+5) -300 ?
<b>Voltage</b>	110V 60Hz /220V 50 Hz
<b>Weight</b>	5.5 kg
<b>Working Size</b>	30x30x12 cm

[www.chinascientific.com](http://www.chinascientific.com)

W Parkway St Denton TX 76201, USA

Email: [info@chinascientific.com](mailto:info@chinascientific.com) | Website: [www.chinascientific.com](http://www.chinascientific.com)