



Portable Sterilizer TRPS-603

www.chinascientific.com | info@chinascientific.com

Description

Portable Sterilizer TRPS-603 is a stainless steel device with an electric/coal heating mode and a Dual scale numerical indication pressure gauge. It is having with a movable interior that is convenient for operation and with a capacity of 18 L and disinfecting Pressure of 0.14-0.165Mpa. It is having with the self-inflating type seal.

Features:

It is made of stainless steel

It has an overpressure auto-discharging system. 0.145 to 0.165MPa

The maximum working temperature is 126 °C

Numeric indication with a dual magnitude Seal that self-inflates

There are two modes: electric heating and electric/coal heating

The movable interior makes operation easier

Specifications :

Capacity	18 L
Dimension	39×39×44 cm
Disinfecting Pressure	0.14-0.165 MPa
Heating Mode	Electric/coal heating
Maximum Working Temperature	126 ?
Power	2 KW
Voltage	110 V 60 Hz /220 V 50 Hz
Weight	12.5 kg

www.chinascientific.com

W Parkway St Denton TX 76201, USA

Email: info@chinascientific.com | Website: www.chinascientific.com