



## Water Distiller TRWD-601

[www.chinascientific.com](http://www.chinascientific.com) | [info@chinascientific.com](mailto:info@chinascientific.com)

## **Description**

Water Distiller TRWD-601 is corrosion and aging resistant and has an intelligent automatic control. It is equipped with a stainless steel bowl with an electric heater. It condenses during the steaming process for the distilled water by electric heating. It has an output water level of 5 liters per hour.

## **Features:**

It is corrosion-resistant and age resistant

Having a smart automatic control

Electric heating causes steam condensation into distilled water

It is outfitted with distillation generating equipment, which includes a condenser and other components

It consists of a stain-free steel electric double distiller

## Specifications :

<b>Output</b>	5 L/ hour
<b>Package</b>	300×300×630 mm
<b>Power</b>	AC220 V/ 4.5 KW
<b>Weight</b>	7 kg

[www.chinascientific.com](http://www.chinascientific.com)

W Parkway St Denton TX 76201, USA  
Email: [info@chinascientific.com](mailto:info@chinascientific.com) | Website: [www.chinascientific.com](http://www.chinascientific.com)